

LEEP PRECISION Integrated System Procedure Preferred Equipment List

DESCRIPTION	PART NUMBER
LEEP PRECISION Interated System with Generator, Smoke Evacuator, and Compact Cart with storage	
<ul style="list-style-type: none"> • Hand or foot activation for operating precision • Elegantly designed for the modern procedure room • Fully integrated LEEP station with high air flow efficiently that captures 99.999% of smoke plumes through triple-stage filtration • Safe and easy to operate 	
LEEP PRECISION Integrated System	
LEEP PRECISION Integrated System 120V	LP-10-120
LEEP PRECISION Integrated System 220V	LP-10-220
LEEP PRECISION Integrated System 220V, China	LP-10-220-CH
LEEP PRECISION Generator	
LEEP PRECISION Generator 120V Standalone	LP-20-120
LEEP PRECISION Generator 220V Standalone	LP-20-220
LEEP PRECISION Generator 220V, Standalone, China	LP-20-220-CH
LEEP PRECISION Smoke Evacuator	
LEEP PRECISION Smoke Evacuator 120V Standalone	LP-30-120
LEEP PRECISION Smoke Evacuator 220V	LP-30-220
LEEP PRECISION Smoke Evacuator 220V, Standalone, China	LP-30-220-CH
LEEP PRECISION Cart	
LEEP PRECISION Cart	LP-40-120
LEEP PRECISION Generator Disposables	
Hand Switch-operated Cautery Pencil with Holster	LP-50-101
Foot Pedal-operated Push Button Cautery with Holster	LP-50-102
Foot Pedal Adapter	LP-50-103
Dual Dispersive Electrodes with Connector	LP-50-201
Dual Dispersive Electrodes without Connector	LP-50-202
Reusable Cable	LP-50-203
LEEP PRECISION Smoke Evacuator Disposables	
Smoke Evacuator Pre-Filter, Single Patient Use (10/bag)	52560
ULPA Filter (1/box)	6082
Smoke Evacuator Reducer, Single Patient Use (10/bag)	6083
Smoke Evacuator Tubing – 6 ft., 1/2 in. outer diameter, Sterile, Single Use (10/box)	6084
Loop, Ball and Cone Electrodes	
Large Radius Loop Electrodes	
2.5 cm W x 1.0 cm D – 12 cm shaft (5/box)	R2510
2.0 cm W x 2.0 cm D – 12 cm shaft (5/box)	R2020
2.0 cm W x 1.5 cm D – 12 cm shaft (5/box)	R2015
2.0 cm W x 1.0 cm D – 12 cm shaft (5/box)	R2010
2.0 cm W x 0.8 cm D – 12 cm shaft (5/box)	R2008
Medium Radius Loop Electrodes	
1.5 cm W x 1.5 cm D – 12 cm shaft (5/box)	R1515
1.5 cm W x 0.7 cm D – 12 cm shaft (5/box)	R1507
1.5 cm W x 0.5 cm D – 12 cm shaft (5/box)	R1505
1.0 cm W x 1.0 cm D – 12 cm shaft (5/box)	R1010
1.0 cm W x 0.7 cm D – 12 cm shaft (5/box)	R1007

Square Electrodes

1.0 cm W x 1.0 cm D – 12 cm shaft (5/box)	S1010
1.0 cm W x 0.7 cm D – 5.5 cm shaft - 15° Bend (5/box)	S1007
1.0 cm W x 0.4 cm D – 5.5 cm shaft - 15° Bend (5/box)	S1004
1.0 cm W x 0.4 cm D – 8 cm shaft (5/box)	V1004

Needle Electrodes

0.8 mm D x 16 mm – 12 cm shaft (5/box)	N1200
0.8 mm D x 16 mm – 5.5 cm shaft (5/box)	N0550
0.25 mm D x 8 mm – 12 cm shaft (5/box)	N2512
0.25 mm D x 5 mm – 5.5 cm shaft (5/box)	N2555

Ball Electrodes

5 mm D – 12 cm shaft (5/box)	B0512
5 mm D – 5.5 cm shaft (5/box)	B0555
3 mm D – 12 cm shaft (5/box)	B0312

Fischer Cone Electrodes

Small 9 mm W x 12 mm D – 15 mm Wire Length	900-150
Medium 12 mm W x 15 mm D – 19 mm Wire Length	900-151
Large 15 mm W x 18 mm D – 23 mm Wire Length	900-152
Medium Extended 12 mm W x 20 mm D – 23 mm Wire Length	900-154
Large Shallow 15 mm W x 8 mm D – 17 mm Wire Length	900-155
Small Wide 9 mm W x 12 mm D – 17.9 mm Wire Length	900-157
Small Wide Extended 9 mm x 20 mm D – 25.8 mm Wire Length	900-158

Insulated Instruments

LEEP Insulated Vaginal Specula

(Note: CooperSurgical has many options for insulated specula – those listed below are the most popular. See the CooperSurgical catalog for more options.)

Graves LEEP Speculum (with smoke tube guide) Large – 125 mm L x 35 mm W	F200
Graves LEEP Speculum (with smoke tube guide) Medium – 113 mm L x 33 mm W	F205
Pederson LEEP Speculum Small – 84 mm L x 16 mm W	F223
Pederson LEEP Speculum Medium – 114 mm L x 22 mm W	F221
Pederson LEEP Speculum Large – 120 mm L x 27 mm W	F220
Vu-More LEEP Speculum Medium – 113 mm L x 33 mm W	F217
Vu-More LEEP Speculum Larger – 125 mm L x 35 mm W	F215

Endocervical Specula

LEEP Endocervical Speculum Regular Length 24 cm – Tip 6 mm	F500
LEEP Endocervical Speculum Medium Length 24 cm – Tip 5 mm	F510

Lateral Vaginal Wall Retractors

LEEP Cer-View Lateral Retractor Medium – 57 mm L x 22 mm W	F410
LEEP Cer-View Lateral Retractor Large – 70 mm L x 22 mm W	F420

Iris Hook

LEEP Iris Hook	F800
----------------	------

Single Tooth Tenaculum

LEEP Single Tooth Tenaculum	F900
-----------------------------	------

Tissue Forceps

LEEP Ring Forceps	F920
LEEP Tissue Forceps	F700

Rulers

LEEP Depth Ruler	F600
LEEP Width Ruler	F610

LEEP Procedure Kits

LEEP Redi-Kit, Single Use, Disposable (5 kits/box)	6061
--	------